- en** Instructions for use - Infusion filter
- de** Gebrauchsanweisung - Infusionsfilter
- cs** Návod k použití - Infuzní filtr
- da** Brugsanvisning - Infusionsfilter
- es** Instrucciones de uso - Filtro de infusión
- fr** Notice d'utilisation - Filtre à perfusion
- it** Istruzioni per l'uso - Filtro di infusione
- nl** Gebruiksaanwijzing - Infusiefilter
- no** Bruksanvisning - Infusjonsfilter
- pl** Instrukcja używania - Filtr infuzyjny
- pt** Instruções de utilização - Filtro de infusão
- sv** Bruksanvisning - Infusionsfilter

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


<https://www.codancompanies.com/download-centre/>


## Infusion filter

suitable for use with compatible pumps as well as gravity feed

**I.V.STAR®** Infusion filter 0.2 micron for infusion therapy for retention of particles and microorganisms as well as elimination of air

**I.V.STAR® PLUS** Infusion filter 0.2 micron (positively charged, endotoxin retentive) for infusion therapy for retention of particles and microorganisms as well as elimination of air

**I.V.STAR® F** Infusion filter 1.2 micron for the administration of fat emulsions for retention of fungi, spores, particles as well as elimination of air

### Colour of housing

**I.V.STAR®** transparent  
**I.V.STAR® PLUS** transparent-green  
**I.V.STAR® F** transparent-blue

### Product characteristics

	1.6 cm <sup>2</sup>	4.5 cm <sup>2</sup>	10 cm <sup>2</sup>
• Filter area :	1.6 cm <sup>2</sup>	4.5 cm <sup>2</sup>	10 cm <sup>2</sup>
• Filling volume :	< 0.4 ml	0.7 ml	2.4 ml
• Flow rate* :			
I.V.STAR®	≥ 2 ml/min	≥ 10 ml/min	≥ 20 ml/min
I.V.STAR® PLUS	≥ 2 ml/min	≥ 20 ml/min	≥ 20 ml/min
I.V.STAR® F	≥ 30 ml/min	≥ 72 ml/min	≥ 100 ml/min
• Pressure resistance :	2 bar [29 psi]	2 bar [29 psi]	2 bar [29 psi]
• Protein binding :	low	low	low

\* Flow rate tested with filtered water at ~23°C / water column 910 mm / with Ø 2.6 mm tubing

### Intended purpose

Filtration for the separation of particles, microorganisms, fungi, spores, endotoxins and air in infusion lines for the transfer of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions by intravenous infusion.

### Indications

Diseases and conditions requiring intravenous infusion of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions.

### Contraindications

The product is not designed and approved for :

- Home use
- Use by unqualified, untrained personnel

Self-application is not permitted and is not considered as intended use.

The product is not intended for the transfer of :

- Blood and blood components
- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions with regard to incompatibilities (interaction material / pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

### Patient target group

According to the intended flow rate and an acceptable dead volume, a corresponding filter surface area must be chosen in accordance with the patient target group.

The available filter size areas are designed as follows :

Filter area 1.6 cm<sup>2</sup> / 4.5 cm<sup>2</sup> : form birth / from birth weight (including premature babies)

Filter area 10 cm<sup>2</sup> : patients of at least 2 years of age and / or a weight of 12 kg

The product can be used for any adult patient and for permitted paediatric patient populations (incl. neonatology) with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

### User

Qualified medical and pharmaceutical personnel only.

For use in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions.

### Clinical benefits

The transfer of fluids, pharmaceuticals and parenteral nutrition solutions by intravenous infusion can bring about therapeutic benefits that improve the health of the patient.

In addition, the use of infusion filters may reduce the risks which may arise from particles, microorganisms, and air.

Infusions of fluids, pharmaceuticals and parenteral nutrition solutions carry implicit risks and should therefore only be carried out if they improve the patient's health.

### General safety instructions

- The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information of this product and compatible medical devices that are used in parallel shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and has to be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids, pharmaceuticals and parenteral nutrition solutions, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material / drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • Prior to the consecutive administration of different drugs with the same infusion filter, it is recommended to perform a benefit-risk assessment of the drug combination, the choice of infusion filter and the planned duration and environmental conditions of use. After each administration of an individual drug, appropriate flushing with neutral, compatible solutions is recommended to avoid potential incompatibilities, maintain the integrity of the infusion filter / line and prevent microbiologic growth. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections of users.

### Notes on safe use of product

- I.V.STAR® / I.V.STAR® PLUS (endotoxin retentive, positively charged) infusion filters are not suitable for administration of blood, blood components, clotting factors, colloidal solutions, suspension, emulsions and undissolved medicaments. • I.V.STAR® F infusion filters are not suitable for administration of blood, blood components, clotting factors, colloidal solutions and undissolved medicaments. • The infusion filters must not be used for sterile filtration. • The product can be used in conjunction with compatible medical devices which are suitable for use with compatible pumps, by gravity or pressure cuff with a maximum generated pressure of up to 50 kPa (0.5 bar). Restrictions to the type of application are possible due to the product structure (e.g. products with an integrated pressure-activated valve cannot be used by gravity feed).

• When used with infusion pumps, inline filters may significantly impact their performance, particularly at low flow rates. They may reduce flow and increase the time required to reach a steady state, especially as particles accumulate on the filter over time. • Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Prior to use, existing closure plugs and disconnectable connections have to be checked for tightness and retightened if necessary. • If equipped with an orange protective cap (vented), replace it with a closure plug (non-vented) before using the medical device. • Infusion filters have to be used as close to the patient as possible to achieve the most effective filtration performance. • Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion container. If higher pressure is applied on the FlowStop protection cap, its function cannot be assured (i.e. protection against leakages). In this case, the FlowStop protection cap has to be replaced by a sterile closure plug immediately after priming. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70% or isopropyl alcohol 70% in combination with chlorhexidine gluconate 2% can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Long-term administration of alcohol- and lipid-based solutions can cause stress cracks. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Blocked infusion filters must not be flushed due to potentially damaging the filter membrane and spreading existing contaminations. Blocked infusion filters must be replaced by a new infusion filter. Blocked infusion filters indicate potential contaminations or incompatibilities of medications. • Avoid excessive injection pressure during injections of solutions into the infusion set/line with integrated infusion filter, especially when using small syringes. Injection pressure that is too high can damage the filter membrane and/or cause leakages in the filter housing. • The opening(s) of the venting filter(s) of the infusion filter should be located in a free and unconcealed position during application. Care must also be taken to ensure that the vent(s) do not come into contact with other substances such as disinfectants, oil or soap from the outside. • In case of administration of light-sensitive pharmaceuticals, an appropriate light-protective cover must be used. Opening(s) of venting filter(s) should nevertheless be located in a free/unconcealed position.

#### Use of product

• Prime the infusion filter with neutral solution in vertical position (for alignment: see exemplary product illustration(s) for reference). Ensure

that all air in the system is eliminated. • Remove remaining air bubbles on the patient side of the infusion filter (serrated side) and/or in the lower connection tube to the patient connection port by tapping/snapping carefully. • Primed infusion filters should be used promptly after priming. Before connection to the patient, check again for potential air bubbles and remove if necessary. • Pharmaceuticals must be completely dissolved to avoid blocked infusion filters.

#### Product compatibility

Compatible with all volumetric infusion and syringe pumps with a pressure alarm limit not exceeding 2.0 bar. See the pump manufacturer's instructions for use.

For single infusion filters with and without integrated tubing connection: All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/EN ISO 80369-7.

#### Residual risks

• Misconnections of the Luer connectors with medical devices of other application areas (e.g. enteral nutrition) are possible. • In the event of handling errors or, in rare cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism and/or under-/overdosing of the infusion. The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

#### Period of use

The international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account.

The product is designed for a period of use of maximum 96 hours.

For non-lipid-containing solutions (use of I.V.STAR®/I.V.STAR® PLUS), a maximum period of use of 96 hours is recommended and 24 hours for lipid-containing solutions (use of I.V.STAR® F).

Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

#### Product disposal

• The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. • Protective medical gloves must be worn during disposal.

#### Notification of serious incidents

• All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay. • Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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