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

Administration

Secondary line

for mixing of cytostatics in particular and administration in conjunction with Cyto-Ad Z[®]/Cyto-Ad-Wing Z[®]

Intended purpose

Transfer of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions by intravenous infusion.

Indications

Diseases and conditions requiring intravenous infusion of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions.

Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained personnel

Self-application is not permitted and is not considered as intended use.

The product is not intended for the transfer of:

- Blood and blood components
- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions with regard to incompatibilities (interaction material/ pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg.

The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

User

Qualified medical and pharmaceutical personnel only. For use in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions.

Clinical benefits

The transfer of fluids, pharmaceuticals and parenteral nutrition solutions by intravenous infusion can bring about therapeutic benefits that improve the health of the patient.

Infusions of fluids, pharmaceuticals and parenteral nutrition solutions carry implicit risks and should therefore only be carried out if they improve the patient's health.

General safety instructions

- The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information of this product and compatible medical devices that are used in parallel shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g., leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and must be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids, pharmaceuticals and parenteral nutrition solutions, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be

observed at all times. • The preparation of pharmaceuticals, especially cytostatics should be carried out under strict and controlled conditions. Appropriate personal protective equipment for the preparation of pharmaceuticals must be used. • The use of protective medical gloves is recommended to prevent any infections of users.

Notes on safe use of product

- The product is intended for gravity application and for infusion by pressure cuff with generated pressures up to 50 kPa (0.5 bar). • The product must be used in conjunction with an infusion set designed for use with compatible pressure infusion pumps, under gravity and if appropriate with a pressure cuff. • For administration, the use in conjunction with CODAN infusion sets Cyto-Ad Z[®]/Cyto-Ad-Wing Z[®] is recommended. • Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Prior to the preparation of pharmaceuticals/ cytostatics, initial priming of the secondary line with neutral solution by gravity feed is recommended. If the secondary line remains unfilled (as air barrier), the drip chamber level drops down after connection of the secondary line and opening of the closure clamp (3/3a/3b). Re-balance the fluid level in the drip chamber of the infusion set (main line) directly after connection of secondary line. • Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion container. If higher pressure is applied on the FlowStop protection cap, its function cannot be assured (i.e. protection against leakages). In this case, the FlowStop protection cap has to be replaced by a sterile closure plug immediately after priming. • For products with a vented closure piercing device (1) (D/E) used with infusion bags, it is recommended to keep the ventilation cap always closed during the preparation of drugs/ cytostatics. • When used with infusion bags, it is recommended to add sterile air to the preparation in order to completely empty the bag and the associated tube section during administration. Alternatively, a product with a vented closure piercing device (1) (D/E) can be used. The ventilation cap of the closure piercing device (1) (D/E) is opened when necessary during administration. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70% or isopropyl alcohol 70% in combination with chlorhexidine gluconate 2% can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Excessive connections of Luer-Lock syringes to the mix and withdrawal adapter (2) can cause damages of the valve. Always connect Luer-Lock syringes to the mix and withdrawal adapter (2) only hand-tight. • Avoid excessive injection pressure during injections of solutions into the infusion device/line, especially when using small syringes. Handling errors can result in leakages. • A fluid filter < 20 µm (according to EN ISO 8536-4/-8) is required for infusion. • Ensure that the infusion containers are positioned at a sufficient height for adequate gravity flow during administration. • Ensure a sufficient height difference between the infusion container with neutral solution (main line) and the infusion container(s) with pharmaceuticals/cytostatics (secondary line(s)) during administration. In order to ensure unobstructed fluid flow, the tubing of the secondary line(s) must not kink or hang down. If necessary, use an infusion bag hanger for the main line. • Long-term administration of alcohol- and lipid-based solutions can cause stress cracks. • After complete administration of pharmaceuticals/cytostatics, close the closure clamp (3/3b) of the secondary line nearby the connection port (7) to the infusion set so that no further air can enter the system.

Use of product

Please also observe the exemplary product illustration(s) (A/B/C/D/E) as well as the graphic brief instructions under notes on product use and all further information/images for product use on and in the product packaging.

Reconstitution of the pharmaceutical/cytostatic und preparation of the secondary line:

1. Check if the sterile packaging is damaged or has been unintentionally opened.

2. Remove the product from the sterile packaging and check for possible damage.
3. In case of infusion bags, secondary lines with non-vented (1) (A/B/C) as well as vented closure-piercing devices (1) (D/E) can be used. If applicable, close the ventilation cap of vented closure-piercing devices (1) (D/E) first. In case of rigid/semi-rigid infusion containers/bottles, a secondary line with a vented closure-piercing device (1) (D/E) is required.
4. Close the closure clamp(s) (3/3a/3b).
5. Remove the protective cap and insert the non-vented closure-piercing device (1) (A/B/C) into the port of the infusion container with a rotating movement or insert the vented closure-piercing (1) (D/E) straight into the uprightly positioned rigid/semi-rigid infusion container/bottle. In case of infusion bags, please ensure that the wall of the bag remains intact during piercing.
- 6a. Initial priming of the secondary line with neutral solution is recommended. In case of infusion bags with sufficient air volume for emptying the bag, the ventilation cap of vented closure-piercing devices (1) (D/E) remains closed. Otherwise open the ventilation cap of vented closure-piercing devices (1) (D/E). Thus, open the closure clamp(s) (3/3a/3b) and prime the secondary line until connection port (7) using gravity. Do not exert any pressure on the infusion container.
- 6b. When used with Cyto-Ad-Wing Z[®] during administration, the secondary line can alternatively remain unfilled (with an air barrier).
7. Close the closure clamp (3/3a) near the closure-piercing device (1) (A/B/C/E) respectively near the mix and withdrawal adapter (2) (D). When using secondary lines with integrated infusion filter (C/E), additionally close the lower closure clamp (3b) near the infusion filter (4). If applicable, reclose the ventilation cap at the vented closure-piercing device (1) (D/E) for a secure transport and storage.
8. Remove the closure plug from the mix and withdrawal adapter (2). Connect a Luer-Lock syringe to the mix and withdrawal adapter (2) and withdraw the required amount of neutral solution for reconstitution of the pharmaceutical/cytostatic (solid) or withdraw a certain amount of neutral solution for discarding, either according to the corresponding amount of pharmaceutical/cytostatic (liquid) to be added or, if necessary, to balance the amount of infusion solution in the infusion container.
9. Reconstitute the pharmaceutical/cytostatic considering the Summary of Product characteristics and the required medication/patient parameters. In order to completely empty the infusion bag and the associated section of tubing during administration, it may be necessary to add sterile air to the preparation, if there is no ventilation at the closure-piercing device (1) (D/E).
10. Connect the Luer-Lock syringe including the pharmaceutical/cytostatic to the mix and withdrawal adapter (2) and insert the solution into the infusion container. The sterile air additionally drawn up remains in the syringe first, if applicable.
11. Draw up again liquid from the infusion container into the syringe to adapt the concentration of the pharmaceutical/cytostatic between syringe and infusion container.
12. Afterwards, insert the complete fluid first and if applicable, followed by the additionally added air into the infusion container.
13. Disconnect the Luer-Lock syringe and close the mix and withdrawal adapter with the present closure plug respectively with a new Connect Z[®]-LOCK closure plug in case of disposal.
14. Label the infusion container with the corresponding medication/patient data. Re-check if the closure clamp(s) (3/3a/3b) and the ventilation cap at the vented closure-piercing device (1) (D/E) is/are correctly closed. A compatible Connect Z[®] Clip can be used to additionally safeguard the connection between the closure piercing device and the infusion container with the pharmaceutical/cytostatic during transport.
15. Place the prepared secondary line with the pharmaceutical/cytostatic into a transport container.

Connection of secondary line(s) to the infusion set CODAN Cyto-Ad Z[®] / Cyto-Ad-Wing Z[®] for administration of pharmaceuticals/cytostatics:

16. Check for possible leakages and if the closure clamp(s) (3/3a/3b) and, if applicable, the ventilation cap of the vented closure-piercing device (1) (D/E) is/are closed.
17. Check the indicated patient/medication information for accuracy.
18. Hang up the infusion container with the pharmaceutical/cytostatic onto the IV pole. Ensure a sufficient height of the infusion

containers and height difference between the infusion container with neutral solution (main line) and the infusion container(s) with the pharmaceuticals/cytostatics (secondary line(s)) for adequate gravity flow during administration.

19. Interrupt the infusion, if already running, at the infusion set (main line) intended for administration by stopping the volumetric infusion pump, if applicable.
20. Close the roller clamp at the infusion set (main line).
21. Close the closure clamp above the drip chamber at the infusion set (main line). If applicable, ensure that the closure clamp(s) of all further unused line(s) is/are closed correctly.
22. Remove the protective cap from the connection port (7) and connect the secondary line (A/B/C/D/E) to the connection port of the corresponding infusion set (main line) intended for administration.
23. When using rigid/semi-rigid infusion/blood containers/bottles, open the ventilation cap of the vented closure-piercing device (1) (D/E), for infusion bags, open as required.
- 24a. Open the closure clamp(s) (3/3a/3b) of the secondary line.
- 24b. If the secondary line is unfilled (with air barrier), the fluid level in the drip chamber of the infusion set (main line) will drop down after connecting the secondary line and opening the closure clamp(s) (3/3a/3b). Re-balance the fluid level in the drip chamber at the infusion set (main line).
25. Open the roller clamp of the infusion set (main line) and set the required medication and patient parameters. If applicable, start the volumetric infusion pump.
26. In case of secondary lines with ONCE-LOCK 2 adapter (6) (B), additionally press the ONCE-LOCK 2 adapter (6) (B) over the connection port (7) to lock the connection.
27. Check the infusion in regular intervals.
28. After complete administration of the infusion container at the secondary line, interrupt the running infusion again by stopping the volumetric infusion pump if applicable and close the roller clamp of the infusion set (main line).
29. Close the closure clamp(s) (3/3b) nearby the connection port (7). The secondary line remains connected to the infusion set (main line).

Product compatibility

The closure-piercing device complies with EN ISO 8536-4 and can be connected with all compatible infusion containers. All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/EN ISO 80369-7.

Residual risks

- Misconnections of the Luer connectors with medical devices of other application areas (e.g., enteral nutrition) are possible.
- In the event of handling errors or, in rare cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism and/or under-/overdosing of the infusion.

The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

Period of use

The international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account.

The product is designed for a period of use of maximum 96 hours. A maximum period of use of 96 hours is recommended for non-lipid-containing solutions and 24 hours for lipid-containing solutions. Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

Product disposal

Do not disconnect any components or medical devices used in parallel. Dispose of the complete system including all connected infusion containers.

The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste.

For disconnection of the main line from the patient connection port and product disposal, appropriate protective measures (e.g. wearing protective medical gloves or personal protective equipment) must be observed.

Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay.

Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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