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

(1) B. Braun Space^{plus} Infusomat[®]/ Infusomat[®] Space

(2) B. Braun Infusomat[®] fmS/fm

Infusion set IVIP® B2 Cyto-Ad Z®/IVIP® B2 Cyto-Ad-Wing Z®

for the administration in particular of cytostatics in conjunction with Connect Z®; equipped with pump segment and Anti-Free-Flow clamp (AFF clamp) suitable for use with the following compatible pumps:

B. Braun Space^{plus} Infusomat®/Infusomat® Space/Infusomat® fmS/Infusomat® fm

Intended purpose

Transfer of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions by intravenous infusion.

Indications

Diseases and conditions requiring intravenous infusion of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions.

Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained staff

Self-application is not permitted and is not considered as intended use.

The product is not intended for the transfer of:

- Blood and blood components
- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions with regard to incompatibilities (interaction material/pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg.

The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

User

Qualified medical and pharmaceutical personnel only.

For use in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions.

Clinical benefits

The transfer of fluids, pharmaceuticals and parenteral nutrition solutions by intravenous infusion can bring about therapeutic benefits that improve the health of the patient.

Infusions of fluids, pharmaceuticals and parenteral nutrition solutions carry implicit risks and should therefore only be carried out if they improve the patient's health.

General safety instructions

• The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g., leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and must be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids, pharmaceuticals and parenteral nutrition solutions, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • In case of consecutive administration of different drugs with the same infusion set (cf. exemplary product illustrations A/B), rinsing with a neutral compatible solution is recommended after each medication to avoid potential incompatibility of drug mixtures, maintain the integrity of the line and prevent microbiological growth. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections of users.

Notes on safe use of product

• Please observe any further instructions for use that are enclosed with this medical device or the instructions for use of any compatible medical devices used in parallel. • The product is designed for use with compatible pressure infusion pumps specified under product compatibility with generated pressures up to 210 kPa (2.1 bar). • The product can be used without using an infusion pump and thus without function of the pump segment as well as the Anti-Free-Flow clamp under gravity and with a pressure cuff with generated pressures of up to 50 kPa (0.5 bar). Restrictions in the type of application are possible due to the product structure (e.g., products with integrated pressure-activated valve cannot be used by gravity feed). • Notes on use with pumps: the user shall ensure that the pump is approved for the administration of fluids, pharmaceuticals and parenteral nutrition solutions and calibrated with the appropriate infusion device before commissioning. Always observe the instructions for use from the pump manufacturer. Deviations from the pump manufacturer's information are possible in these instructions for use and shall be observed. Only open the roller clamp during pump operation! • Only insert the Anti-Free-Flow clamp into the B. Braun Space^{plus} Infusomat®/Infusomat® Space when the opening lever is

open. • The Anti-Free-Flow clamp remains unused for B. Braun Infusomat® fmS/fm pumps. • The use of CODAN secondary lines Connect Z® and Koppelsystem in conjunction with infusion sets IVIP® B2 Cyto-Ad Z®/IVIP® B2 Cyto-Ad-Wing Z® is recommended. • Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Prior to use, existing closure plugs and disconnectable connections must be checked for tightness and retightened if necessary. • If equipped with an orange protective cap (vented), replace it with a closure plug (non-vented) before using the medical device. • A fluid filter < 20 µm (according to EN ISO 8536-8) is required for infusion. • Prior to the administration of pharmaceuticals/cytostatics, initial priming of the product with neutral solution by gravity feed is recommended. • Do not stretch/extend the pump segment (6a-c) excessively. • Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion container. If higher pressure is applied on the FlowStop protection cap, its function cannot be assured (i.e., protection against leakages). In this case, the FlowStop protection cap has to be replaced by a sterile closure plug immediately after priming. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70 % or isopropyl alcohol 70 % in combination with chlorhexidine gluconate 2 % can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Ensure a sufficient height difference between the infusion container with neutral solution (main line) and the infusion container(s) with pharmaceuticals/cytostatics (secondary line(s)) during administration. In order to ensure unobstructed fluid flow, the tubing of the secondary line(s) must not kink or hang down. If necessary, use an infusion bag hanger for the main line. • In case of consecutive administration of different pharmaceuticals/cytostatics, rinsing with sufficient neutral solution (min. 50 ml in between of two different pharmaceuticals/cytostatics and min. 100 ml after the last administration of pharmaceutical/cytostatic) is recommended. • Long-term administration of alcohol- and lipid-based solutions can cause stress cracks. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Avoid excessive injection pressure during injections of solutions into the infusion device/line, especially when using small syringes. Handling errors can result in leakages. • Before removing the pump segment (6a-c) from the infusion pump, stop the running infusion, close the roller clamp (9) and disconnect the patient connection port (11) from the patient access to avoid any free-flow. • LIGHT-PROTECTED product: The protective cover of the drip chamber can be pulled down to check the infusion or, if necessary, to connect the drop detector of the volumetric infusion pump.

Use of product

Please also observe the exemplary product illustration(s) (A/B) as well as the graphic brief instructions under notes on product use and further information/images for product use on and in the product packaging.

Preparation of the infusion set:

1. Check if the sterile packaging is damaged or unintentionally opened.
 2. Remove the product from the sterile packaging and check for possible damage.
 3. Close the closure clamp (2) above the drip chamber (5a/b).
 4. Tighten all closure plugs on the connection ports (3), so they are hand-tight.
 5. Close the roller clamp (9).
 6. In case of rigid/semi-rigid infusion containers/bottles, an infusion device with a vented closure-piercing device (1) (B) is required. In case of infusion bags, infusion devices with non-vented (1) (A) as well as vented closure-piercing devices (1) (B) can be used. Close the ventilation cap of vented closure-piercing devices (1) (B) first.
 7. Remove the protective cap and insert the vented closure-piercing (1) (B) straight into the upright-positioned rigid/semi-rigid infusion container/bottle or insert the non-vented closure-piercing device (1) (A) into the port of the infusion bag with a rotating movement. Please ensure that the wall of the bag remains intact during piercing.
 8. Hang up the infusion container onto the IV pole.
 9. In the case of rigid/semi-rigid infusion containers/bottles, open the ventilation cap of the vented closure-piercing device (1) (B). For infusion bags, the ventilation cap of the vented closure-piercing device (1) (B) remains closed.
 - 10a. While holding the drip chamber (5a) (A) at an angle upwards, first open the closure clamp (2) and then the roller clamp (9) to fill the drip chamber from top to bottom to approximately 1 cm below the level ring. Once the level has been reached, close the roller clamp (9) and return the drip chamber (5a) back into its vertical starting position.
 - 10b. Regarding products equipped with Wing-Valve drip chamber (5b) (B), open the closure clamp (2) and fill the drip chamber (5b) to approximately 1 cm below the level ring (4) or regulate the fluid level by pressing the Wing-Valve wings together while the drip chamber (5b) remains in its vertical position all the time.
 11. Open the roller clamp (9) and prime the remaining infusion set up to the patient connection port (11) so that the product is applied without any air inclusions. To avoid any air inclusions, it is recommended to hold the pump segment (6a-c) upside down and to fill it without stretching the pump segment from above.
 12. Close the roller clamp (9).
 13. In case of gravity application: proceed according to point 22.
When using with pump: follow the pump manufacturer's instructions for use.
 - 13a. Ensure that a maximum distance between the infusion container and infusion pump is established.
 - 13b. Open the door of the infusion pump.
- Insert the pump segment and the Anti-Free-Flow clamp (AFF clamp) into
B. Braun Space^{plus} Infusomat®/Infusomat® Space (cf. Notes on product use fig. (1)):
 14. Insert the product with pump segment (6a-c) straight from right to left.

15. Connect the upstream connector (6a) of the pump segment to the fixation device on the right side of the delivery element of the pump. The arrow on the upstream connector (6a) must be visible. Ensure that the tubing (8) in front of the upstream connector (6a) is inserted into the tubing guidance channel of the infusion pump.
16. Place the silicone tubing (6b) over the delivery element of the pump without excessive elongation.
17. Then, connect the downstream connector (6c) of the pump segment to the fixation device on the left side of the delivery element. Make sure that the silicone tubing (6b) with CODAN lettering is readable/aligned/not twisted.
18. Insert the AFF clamp (7) into the opening at the infusion pump until the opening lever locks in. The arrows on the upper side of the AFF clamp indicate the direction of insertion. The warning light on the pump goes out.
19. After that, insert the tubing (8) into the tubing guidance channel of the air and pressure sensors of the pump.
20. Finally, check if the tubing (8) is properly inserted into the tubing guidance channels of the infusion pump. Close the door of the infusion pump while making sure that no tubing (8) is compressed by the infusion pump door. The roller clamp (9) must be positioned on the patient connection port (11) site.
21. If applicable, attach the drip detector of the pump to the ring of the drip chamber (5a/b).

Insert the pump segment into **B. Braun Infusomat® fmS/fm** (cf. Notes on product use fig. (2)):

14. Insert the product with pump segment (6a-c) straight from top to bottom. First, insert the tubing (8) into the tubing guidance device at the top of the pump handle.
15. Connect the upstream connector (6a) of the pump segment to the fixation device on the top of the delivery element of the pump. The arrow on the upstream connector (6a) must be visible.
16. Place the silicone tubing (6b) over the delivery element of the pump without excessive elongation. Make sure that the silicone tubing (6b) is not yet squeezed into the constriction (anti-free-flow protection) directly below the delivery element.
17. Then, connect the downstream connector (6c) of the pump segment to the fixation device on the bottom of the delivery element. Make sure that the silicone tubing (6b) with CODAN lettering is readable aligned from bottom to top/not twisted.
18. Insert the tubing (8) below the pump segment (6a-c) into the tubing guidance channel of the air sensor at the lower part of the infusion pump.
19. The AFF clamp (7) remains below the infusion pump without being used. Ensure that the AFF clamp (7) remains open.
20. Finally, check if the tubing (8) is properly inserted into the tubing guidance channels of the infusion pump. Close the door of the infusion pump and make sure that no tubing (8) is compressed by the infusion pump door. The roller clamp (9) must be positioned on the patient connection port (11) site.
21. If applicable, attach the drip detector of the pump to the ring of the drip chamber (5a/b).

Start the infusion:

22. Remove the protective cap from the patient connection port (11) and connect the patient connection port (11) to the patient access.
23. Adjust the appropriate required medication and patient parameters, open the roller clamp (9), and start the infusion.
24. Check the infusion at regular intervals.

Connection of secondary line(s) to the infusion set for the administration of pharmaceuticals/cytostatics:

25. After checking for intactness, remove the prepared secondary line with the connected infusion container including the pharmaceutical/cytostatic from the transport container. Check if the closure clamp(s) is/are closed.
26. Check the indicated patient/medication information for accuracy.
27. Hang up the infusion container with the pharmaceutical/cytostatic onto the IV pole.
28. Interrupt the infusion, if in progress, by stopping the volumetric infusion pump, if applicable.
29. Close the roller clamp (9).
30. Close the closure clamp (2) above the drip chamber (5a/b) of the infusion set. If applicable, ensure that the closure clamp(s) of all further unused line(s) is/are closed correctly.
- 31a. Remove the closure plug from one connection port (3) of the infusion set and connect the secondary line with infusion container. Open the closure clamp(s) of the secondary line.
- 31b. If the secondary line remained unfilled (with air barrier), the drip chamber level will drop down after connection of the secondary line and opening of the closure clamp. The fluid level in the drip chamber (5a/b) has to be rebalanced.
32. Open the roller clamp (9) and set the required medication and patient parameters. If applicable, start the volumetric infusion pump. Ensure a sufficient height difference between the infusion container with neutral solution (main line) and the infusion container(s) with the pharmaceutical/cytostatic (secondary line).
33. Check the infusion at regular intervals.
34. After complete administration of the pharmaceutical/cytostatic, interrupt the infusion again by stopping the volumetric infusion pump, if applicable, and close the roller clamp (9).
35. Close the closure clamp(s) of the drained secondary line near the connection port (3).

36. Open the closure clamp (2) above the drip chamber (5a/b). If the drip chamber level drops down, refill the drip chamber (5a/b) up to the level ring (4). Open the roller clamp (9) and start the volumetric infusion pump, if applicable, to rinse with sufficient neutral solution (min. 50 ml in between of two different pharmaceuticals/cytostatics and min. 100 ml after the last administration of the pharmaceutical/cytostatic).
37. In case of consecutive administration of different pharmaceuticals/cytostatics via secondary lines, repeat steps 25. to 36.. Ensure that the closure clamp(s) of all unused secondary line(s) as well as the closure clamp (2) of the main line are closed.

Completion of the infusion:

38. After completion of the infusion and the final flushing of the main line, close all closure clamp(s) and the roller clamp (9) on the entire infusion system.
39. First, disconnect the patient connection port (11) from the patient access. Do not disconnect any further components or medical devices used in parallel.
40. If applicable, open the door of the infusion pump.
41. Additionally for B. Braun Space^{plus} Infusomat[®]/Infusomat[®] Space: Remove the AFF clamp (7) by pressing the opening lever downwards until it locks into place and the AFF clamp (7) is released.
42. Unlock the connectors (6a/6c) from the fixation devices and remove the product with the pump segment (6a-c) from the infusion pump.
43. Dispose of the complete system including all connected secondary line(s) with infusion container(s).

Product compatibility

The product is compatible with the following volumetric infusion pumps:

B. Braun Space^{plus} Infusomat[®]/Infusomat[®] Space/Infusomat[®] fmS/Infusomat[®] fm

Specification of occlusion pressure alarm* with CODAN products:

Pump	Pressure level	CODAN PVC (TOTM) tubing		CODAN PVC-free tubing (green line [®])			
		Time to occlusion pressure alarm at rate [mm:ss]		Max. bolus volume [ml]	Time to occlusion pressure alarm at rate [mm:ss]		Max. bolus volume [ml]
		75 ml/h	200 ml/h		75 ml/h	200 ml/h	
Space ^{plus}	Level 1	0:06	0:02	0.1	0:07	0:03	0.2
Infusomat [®]	Level 9	0:27	0:10	0.5	0:21	0:09	0.3
Infusomat [®]	Level 1	0:10	0:04	0.0	0:06	0:03	0.0
Space	Level 9	0:34	0:12	0.1	0:21	0:07	0.0
Infusomat [®]	medium	0:16	0:06	0.3	0:13	0:05	0.3
fmS	high	0:23	0:08	0.5	0:18	0:07	0.4
Infusomat [®]	medium	0:07	0:03	0.2	0:06	0:03	0.2
fm	high	0:18	0:07	0.4	0:12	0:05	0.3

* Testing in accordance to IEC 60601-2-24. The type of the infusion line, any integrated components (e.g., infusion filter), temperature, as well as the viscosity of the infusion fluid also influence the time until occlusion pressure alarm and the maximal bolus volume next to the pressure level.

The closure-piercing device complies with EN ISO 8536-8 and can be connected with all compatible infusion containers.

All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/EN ISO 80369-7.

Residual risks

- Misconnections of the Luer connectors with medical devices of other application areas (e.g., enteral nutrition) are possible.
- In the event of handling errors or, in rare cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism and/or under-/overdosing of the infusion.

The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

Period of use

The international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account.

The product is designed for a period of use of maximum 96 hours.

A maximum period of use of 96 hours is recommended for non-lipid-containing solutions and 24 hours for lipid-containing solutions.

Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

Product disposal

Do not disconnect any components or medical devices used in parallel. Dispose of the complete system including all connected infusion containers.

The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste.

During disconnection from the patient connection port and product disposal, appropriate protective measures (e.g., wearing of protective medical gloves or personal protective equipment) must be observed.

Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay. Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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