en Instructions for use

Closed Male Luer adapter
CLOSE-LOCK

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


1b

2

Closed Male Luer adapter CLOSE-LOCK

Intended purpose

Closure of a Luer-Lock connector as well as injection, aspiration and transfer via septum into infusion/transfusion lines for the transfer of fluids, pharmaceuticals, parenteral nutrition solutions or blood and blood components by intravenous infusion.

Indications

Diseases and conditions requiring intravenous infusion of fluids, pharmaceuticals and parenteral nutrition solutions or blood and blood components.

Contraindications

The product is not designed and approved for:

- Home use
- Use by unqualified, untrained personnel

Self-application is not permitted and is not considered as intended use.

The product is not intended for the transfer of:

- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions or blood and blood components with regard to incompatibilities (interaction material/pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

Clinical benefits

The transfer of fluids, pharmaceuticals and parenteral nutrition solutions or blood and blood components by intravenous infusion can bring about therapeutic benefits that improve the health of the patient.

Infusions of fluids, pharmaceuticals and parenteral nutrition solutions and transfusions of blood and blood components carry implicit risks and should therefore only be carried out if they improve the patient's health.

General safety instructions

• The product is a medical device. • The product shall only be used by qualified medical and pharmaceutical personnel in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage) or performance changes, the product shall not be used on the patient and must be replaced. The defective product is to be provided to the manufacturer in compliance with

measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids, pharmaceuticals and parenteral nutrition solutions or blood and blood components, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions or blood and blood components shall be taken into account. In particular, blood group compatibility between donor and recipient as well as the EN ISO 1135-5 recommendation not to transfer platelet components under pressure shall be considered. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • The function of the product shall be checked by qualified medical/pharmaceutical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections of users.

Notes on safe use of product

• Please observe any further instructions for use that are enclosed with this medical device or the instructions for use of any compatible medical devices used in parallel. • The product can be used in conjunction with compatible medical devices which are suitable for use with compatible pumps (max. 200 kPa/2.0 bar), by gravity or pressure cuff with a maximum generated pressure of up to 50 kPa (0.5 bar). • Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • A protective cap is provided with the product in the packaging. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70% or isopropyl alcohol 70% in combination with chlorhexidine gluconate 2% can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • The product is designed for use with Luer-Lock connectors/syringes and can be used for the preparation and administration of pharmaceuticals. • Do not use cannulas or similar with the product. • When using the product at the distal end of a secondary line Connect Z[®]/Koppelsystem equipped with a ONCE-LOCK 2 Male Luer-Lock connector, the locked connection ONCE-LOCK 2 to the female Luer-Lock connector of Cyto-Ad Z[®]/Cyto-Ad-Wing Z[®]

is no longer possible. • For connection/disconnection, the product is held on the side facing the user (gripping surface of the Male Luer-Lock connector). • The Female Luer-Lock connector cannot be disconnected once it has been connected to a Male Luer-Lock connector/syringe. • A simultaneous push and rotary movement is required to connect the Male Luer-Lock connector. The septum of the valve in the Male Luer of the product opens automatically when a Female Luer-Lock connector is connected, so that fluid can flow in both directions. When the Female Luer-Lock connector is disconnected, the septum of the Male Luer valve closes automatically. • The product can be primed by using a Female Luer-Lock connector. The priming volume is 0.1 ml. • The white protective cap does not activate the Male Luer-Lock connector of the product. • Connect the white protective cap to the Male Luer-Lock connector for protection during transport. • Long-term administration of alcohol- and lipid-based solutions can cause stress cracks. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Avoid excessive injection pressure during injections of solutions into the infusion/transfusion device/line, especially when using small syringes. Handling errors can result in leakages. • When administering light-sensitive pharmaceuticals, the product shall also be protected against light exposure by suitable measures.

Product compatibility

All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/EN ISO 80369-7.

Residual risks

The following residual risks are known when using the product on patients:

- Misconnections of the Luer connectors with medical devices of other application areas (e.g., enteral nutrition) are possible. • In the event of handling errors or, in rare

cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism and/or under-/overdosing of the infusion. • The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

Period of use

The international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions or blood and blood components shall be taken into account.

The product is designed for a period of use of maximum 96 hours.

A maximum period of use of 96 hours is recommended for non-lipid-containing solutions and 24 hours for lipid-containing solutions.

A maximum period of use of 12 hours should not be exceeded for blood and blood components.

Taking into account the individual benefit-risk ratio of each patient and the summary of product characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

Product disposal

- The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. • Protective medical gloves must be worn during disposal.

Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay. Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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