
- en** **Instructions for use - Transfusion set**  
for gravity application
- de** **Gebrauchsanweisung - Transfusionsgerät**  
für die Verwendung unter Schwerkraft
- cs** **Návod k použití - Transfuzní souprava**  
pro gravitační aplikaci
- da** **Brugsanvisning - Transfusions sæt**  
til anvendelse under tyngdekraft
- es** **Instrucciones de uso - Equipo de transfusión**  
para aplicación por gravedad
- fr** **Notice d'utilisation - Nécessaire pour transfusion**  
pour application par gravité
- it** **Istruzioni per l'uso - Set per trasfusione**  
per uso tramite sfruttamento della forza di gravità
- nl** **Gebruiksaanwijzing - Transfusiesysteem**  
voor gebruik met zwaartekracht
- no** **Brukerveiledning - Transfusjonssett**  
for bruk under falltilførsel
- pl** **Instrukcja używania - Zestaw do transfuzji**  
do podawania grawitacyjnego
- pt** **Instruções de utilização - Sistema de transfusão**  
para aplicação sob gravidade
- sv** **Bruksanvisning - Transfusionsaggregat**  
för användning med tyngdkraft

en Valid for  de Gültig für  cs Platné pro  da Gyldigt til  es Válida para  fr Valable pour  it Valido per  nl Geldig voor  no Gyldig for  
 pl Obowiązuje dla  pt Valido para  sv Gäller för

45.5560 B88

(en) Exemplary product illustration(s) - deviations in product configuration are possible! (de) Exemplarische Produktabbildung(en) - Abweichungen in der Produktkonfiguration sind möglich! (cs) Příkladová zobrazení produktu - konfigurace produktu se může lišit! (da) Produktillustrationseksempler - afvigelser i produktkonfiguration kan forekomme! (es) Es posible que haya diferencias entre la ilustración o las ilustraciones de ejemplo del producto y la configuración del producto. (fr) Illustration(s) du dispositif à titre d'exemple - Des différences de configuration du dispositif sont possibles. (it) Illustrazione/i del prodotto a titolo esemplificativo; è possibile che vi siano differenze nella configurazione del prodotto! (nl) Voorbeeldillustratie(s) van product - de configuratie van uw product kan afwijken! (no) Produktillustrasjon(er) - forskjellige produktkonfigurasjoner kan forekomme. (pl) Przykładowa(-e) ilustracja(-e) produktu - zastrzega się możliwość zmian w konfiguracji produktu! (pt) Ilustração(ões) exemplificativa(s) do(s) produto(s) - são possíveis desvios na configuração do produto! (sv) Exempelbild(-er) av produkten - avvikelser i produktens konfiguration är möjliga!

- (en) Closure-piercing device with protective cap (de) Einsteckteil mit Schutzkappe (cs) Uzavírací-propichovací díl s ochrannou krytkou (da) Spike-enhed med beskyttelseshætte (es) Dispositivo de perforación de cierre con caperuzas de protección (fr) Dispositif de perforation du bouchon avec capuchon de protection (it) Dispositivo di perforazione della chiusura con cappuccio protettivo (nl) Afsluitprikker met beschermkapje (no) Gjennomhullingsenhet med beskyttelseshette (pl) Łącznik kolcowy z kapturkiem ochronnym (pt) Dispositivo de perfuração com tampa de proteção (sv) Insticksdel med skyddshätta
- (en) Level ring (de) Niveauring (cs) Hladinový kroužek (da) Niveauring (es) Anillo de nivel (fr) Repère circulaire de niveau (it) Anello di livello (nl) Ring voor peil (no) Nivåmerke (pl) Pierścień poziomy (pt) Anel de nível (sv) Nivåring
- (en) Drip chamber (de) Tropfkammer (cs) Kapací komůrka (da) Dråbekammer (es) Cámara de goteo (fr) Chambre compte-gouttes (it) Camera di gocciolamento (nl) Druppelkamer (no) Dryppkammer (pl) Komora kroplowa (pt) Câmara de gotejamento (sv) Droppkammare
- (en) Roller clamp (de) Rollenklemme (cs) Tlačka (da) Rulleklemme (es) Abrazadera de rodillos (fr) Clamp à roulette (it) Roller clamp (nl) Rolklem (no) Rulleklemme (pl) Zacisk rolkowy (pt) Grampo de rolos (sv) Rullklämma
- (en) Tubing (de) Schlauch (cs) Hadička (da) Slangar (es) Tubo flexible (fr) Tubulure (it) Tubazione (nl) Slangen (no) Slange (pl) Przewód (pt) Tubos (sv) Slangar
- (en) Patient connection port (de) Patientenanschluss (cs) Připojovací port pacienta (da) Patienttilslutningsport (es) Puerto de conexión del paciente (fr) Port de connexion patient (it) Porta di connessione del paziente (nl) Patiëntaansluitpoort (no) Pasientforbindelsespport (pl) Port przyłączeniowy pacjenta (pt) Porta de conexão do doente (sv) Patientanslutningsport

15°C 25°C

30% 60%

## Transfusion set

for gravity application

### Intended purpose

Transfer of blood and blood components by intravenous infusion. Isotonic electrolyte solution can additionally be used for priming and rinsing the product.

### Indications

Diseases and conditions requiring intravenous infusion of blood and blood components.

### Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained staff

Self-application is not permitted and is not considered intended use.

The product is not intended for the transfer of:

- Liquids, pharmaceuticals or parenteral nutrition solutions except those specified under intended use
- Enteral nutrition solutions

The product is not intended for use with pressure infusion pumps. The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of blood/blood components/electrolyte solutions with regard to incompatibilities (interaction material/ pharmaceutical and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention has to be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

### Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg. The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

### User

Qualified medical personnel only. For use in hospitals and specialised practices with conventional environmental and hygiene conditions.

### Clinical benefits

Life-saving and/or health-improving therapeutic benefits to patients. The transfusion of blood and blood components carries implicit risks and should therefore only be carried out if it improves patient's health.

### General safety instructions

- The product is a medical device. • The instructions for use shall be observed before using the product. • If the sterile packaging is damaged or unintentionally opened before use, the product has to be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and has to be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of blood and blood components, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of blood and blood components shall be taken into account. In particular, blood group compatibility between donor and recipient as well as the EN ISO 1135-5 recommendation not to transfer platelet components under pressure shall be considered. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities and the current Summary of Product Characteristics. • Generally accepted hygiene measures and working methods have to be observed at all times. • Use of protective medical gloves is recommended to prevent any infections of users.

### Notes on safe use of product

- The product is intended for gravity application and for infusion by pressure cuff with generated pressures up to 50 kPa (0.5 bar) only.
- Before using the product, compare the permissible fluid intake of the

patient with the suitability of the product. • A complete wetting of the fluid filter in the drip chamber with isotonic electrolyte solution is recommended prior to the administration of blood and blood components. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Existing closure plugs have to be checked on tightness and retightened if necessary. • Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion/blood container. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70 % or isopropyl alcohol 70 % in combination with chlorhexidine gluconate 2 % can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Excessive forces on Luer-Lock connections can lead to stress cracks. • If equipped with an orange protective cap (vented), replace it with a closure plug (non-vented).

### Use of product

Please also observe the exemplary product illustration(s) as well as further information/images for product use on the product packaging.

General use (cf. exemplary product illustration(s)):

1. Check if the sterile packaging is damaged or unintentionally opened.
2. Remove the product from the sterile packaging and check for possible damage.
3. Close the roller clamp (4).
4. Remove the protective cap and insert the non-vented closure-piercing device (1) into the infusion/blood bag with a rotating movement. Please ensure that the wall of the blood bag remains intact during piercing.
5. In case of rigid/semi-rigid infusion or blood containers/bottles, a transfusion device with a vented closure-piercing device is required. Ensure that the ventilation cap on the closure-piercing device (1) is closed first. Insert the closure-piercing device (1) straight into the uprightly positioned rigid/semi-rigid infusion or blood container/ bottle.
6. Fill the drip chamber (3) up to the level ring (2).
7. In the case of rigid/semi-rigid infusion or blood containers/bottles, open the ventilation cap on the closure-piercing device (1).
8. Open the roller clamp (4) and prime the remaining transfusion device up to the patient connection port (6) so that the transfusion device is applied without any air inclusions.
9. Close the roller clamp (4) and connect the patient connection port (6) to the patient access.
10. For starting the transfusion, open the roller clamp (4) and adjust the required medication and patient parameters.
11. Check the transfusion at regular intervals.

### Product compatibility

The closure-piercing device complies with EN ISO 1135-4 and can be connected with all compatible infusion/blood containers. All Luer and Luer-Lock connectors are compatible according to EN 1707 and EN 20594-1.

### Residual risks

Due to leakage or malfunction, biological, physical or chemical contaminations, non-sterility, misconnections, degradations of or damage to blood components, changed volume flow, chemical or physical interactions, the following possible health risks can occur:

- Infections
- Skin irritations/skin injuries
- Air embolism
- General degradation of health status
- Toxic effects

The patient or their legally responsible person should be informed about potential health risks prior to treatment.

### Period of use

The international and national guidelines for the transfer of blood and blood components shall be taken into account. A maximum period of use of 12 hours should not be exceeded. Taking into account the individual benefit-risk ratio of each patient, a longer period of use can be considered by the physician.

## Product disposal

The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. Insert the closure-piercing device into the notch on the roller clamp to avoid any injuries. Protective medical gloves have to be worn during disposal.

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