en Instructions for use - Transfusion set IVIP® 708/718

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
MD

Transfusion set IVIP® 708/718

with Anti-Free-Flow clamp (AFF clamp) suitable for use with the following compatible pumps
 CODAN A718V Plus
 CODAN A718V
 CODAN ARGUS 718 V
 CODAN ARGUS 708 V

Intended purpose

Transfer of blood and blood components by intravenous infusion. Isotonic electrolyte solution can additionally be used for priming and rinsing the product.

Indications

Diseases and conditions requiring intravenous infusion of blood and blood components.

Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained staff

Self-application is not permitted and is not considered intended use.

The product is not intended for the transfer of:

- Liquids, pharmaceuticals or parenteral nutrition solutions except those specified under intended use
- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of blood/blood components/electrolyte solutions with regard to incompatibilities (interaction material/ pharmaceutical and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention has to be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg. The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

User

Qualified medical personnel only. For use in hospitals and specialised practices with conventional environmental and hygiene conditions.

Clinical benefits

The transfer of blood and blood components by intravenous infusion can bring about therapeutic benefits that improve the health of the patient. Transfusions of blood and blood components carry implicit risks and should therefore only be carried out if they improve the patient's health.

General safety instructions

• The product is a medical device. • The instructions for use shall be observed before using the product. If applicable further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product has to be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and has to be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of blood and blood components, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of blood and blood components shall be taken into account. In particular, blood group compatibility between donor and recipient as well as the EN ISO 1135-5 recommendation not to transfer platelet components under pressure shall be considered. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on

request for all medical devices manufactured by CODAN. • Generally accepted hygiene measures and working methods have to be observed at all times. • Use of protective medical gloves is recommended to prevent any infections of users.

Notes on safe use of product

• The product is designed for use with compatible pressure infusion pumps with generated pressures up to 200 kPa (2 bar). • The product can also be used without using the Anti-Free-Flow clamp under gravity and with a pressure cuff with generated pressures of up to 50 kPa (0.5 bar). Restrictions in the type of application are possible due to the product structure (e.g. products with integrated pressure-activated valve cannot be used by gravity feed). • Notes on use with pumps: the user shall ensure that the pump is approved for the administration of blood and blood components and calibrated with the appropriate transfusion device before commissioning. Always observe the instructions for use from the pump manufacturer. Only open the roller clamp during pump operation! • Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Prior to use, existing closure plugs and disconnectable connections have to be checked for tightness and retightened if necessary. • If equipped with an orange protective cap (vented), replace it with a closure plug (non-vented) before using the medical device. • A complete wetting of the fluid filter in the drip chamber with isotonic electrolyte solution is recommended prior to the administration of blood and blood components. • Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion container. If higher pressure is applied on the FlowStop protection cap, its function cannot be assured (i.e. protection against leakages). In this case, the FlowStop protection cap has to be replaced by a sterile closure plug immediately after priming. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70 % or isopropyl alcohol 70 % in combination with chlorhexidine gluconate 2 % can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Avoid excessive injection pressure during injections of solutions into the transfusion device/line, especially when using small syringes. Handling errors can result in leakages.

Use of product

Please also observe the exemplary product illustration(s) as well as the graphic brief instructions under notes on product use and further information/images for product use on and in the product packaging.

General use (cf. exemplary product illustrations A/B):

1. Check if the sterile packaging is damaged or unintentionally opened.
2. Remove the product from the sterile packaging and check for possible damage.
3. Close the roller clamp (7).
4. Remove the protective cap and insert the non-vented closure-piercing device (1) into the infusion/blood bag with a rotating movement. Please ensure that the wall of the infusion/blood bag remains intact during piercing.
5. In case of rigid/semi-rigid infusion or blood containers/bottles, a transfusion device with a vented closure-piercing device is required. Ensure that the ventilation cap on the closure-piercing device (1) is closed first. Insert the closure-piercing device (1) straight into the uprightly positioned rigid/semi-rigid infusion or blood container/bottle.
6. Fill the drip chamber (5) up to the level ring (4).
7. In the case of rigid/semi-rigid infusion or blood containers/bottles, open the ventilation cap on the closure-piercing device (1).
8. Open the roller clamp (7) and prime the remaining transfusion set up to the patient connection port (9) so that the transfusion set is applied without any air inclusions. Close the roller clamp (7).
9. Open the door of the infusion pump.
10. Insert the tubing into the infusion pump (a). Roller clamp (7) has to be positioned on the patient connection port (9) site.

11. Insert the Anti-Free-Flow clamp (6) into the access port (b). Close the door of the infusion pump (c).
12. Connect the patient connection port (9) to the patient access.
13. For starting the transfusion, open the roller clamp (7) and adjust the infusion pump according to required medication and patient parameters.
14. Check the transfusion at regular intervals.

Use of product configurations with connection port(s) above the drip chamber (exemplary product illustration B): A secondary line can be connected to the connection port (3) above the drip chamber (5) for successive administration of several blood containers or successive administration of blood/blood components and rinsing solution. The respective unused line(s) shall be closed.

Product compatibility

The product is compatible with the following volumetric infusion pumps:
 CODAN A718V Plus
 CODAN A718V
 CODAN ARGUS 718 V
 CODAN ARGUS 708 V

The closure-piercing device complies with EN ISO 1135-5 and can be connected with all compatible infusion/blood containers. All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/ EN ISO 80369-7.

Residual risks

• Misconnections of the Luer connectors with medical devices of other application areas (e.g. enteral nutrition) are possible. • In the event of handling errors or, in rare cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism. The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

Period of use

The international and national guidelines for the transfer of blood and blood components shall be taken into account. A maximum period of use of 12 hours should not be exceeded. Taking into account the individual benefit-risk ratio of each patient, a longer period of use can be considered by the physician.

Product disposal

The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. The closure-piercing device can be inserted into the notch on the roller clamp to avoid any injuries. Protective medical gloves have to be worn during disposal.

Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay. Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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