- en** **Instructions for use - Irrigation set**  
for urological and gynaecological applications
- de** **Gebrauchsanweisung - Spülset**  
für urologische und gynäkologische Anwendungen
- cs** **Návod k použití - Irigační sada**  
pro urologická a gynekologická použití
- da** **Brugsanvisning - Irrigationssæt**  
til urologisk og gynækologisk anvendelse
- es** **Instrucciones de uso - Set de irrigación**  
para aplicaciones urológicas y ginecológicas
- fr** **Notice d'utilisation - Set d'irrigation**  
pour applications urologiques et gynécologiques
- it** **Istruzioni per l'uso - Set di irrigazione**  
per applicazioni urologiche e ginecologiche
- nl** **Gebruiksaanwijzing - Irrigatieset**  
voor urologische en gynaecologische toepassingen
- no** **Brukerveiledning - Irrigasjonssett**  
til urologiske og gynekologiske behandlinger
- pl** **Instrukcja używania - Zestaw irygacyjny**  
do zastosowań urologicznych i ginekologicznych
- pt** **Instruções de utilização - Sistema de irrigação**  
para aplicações urológicas e ginecológicas
- sv** **Bruksanvisning - Irrigationsaggregat**  
för urologiska och gynekologiska användningar

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
## Irrigation set

for urological and gynaecological applications

### Intended purpose

Transfer of fluids and pharmaceuticals for irrigation during urological and gynaecological applications

### Indications

Therapeutic or diagnostic urological applications in the urinary bladder and prostate, as well as gynaecological applications in the uterus that require the transfer of fluids and pharmaceuticals.

### Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained staff
- Self-application is not permitted and is not considered as intended use.

The current Summary of Product Characteristics of the pharmaceutical manufacturers with regard to incompatibilities (interaction material/ pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

### Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg.

The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of urethral/uterine anatomy, appropriateness of the solution being used, and the duration of use.

### User

Qualified medical personnel only. For use in hospitals and specialised practices with conventional environmental and hygiene conditions.

### Clinical benefits

The transfer of fluids and pharmaceutical solutions for irrigation in urological and gynaecological applications can bring about clinical benefits that improve the health of the patient.

Irrigation in urological and gynaecological applications carries implicit risks and should therefore only be carried out if it improves the patient's health.

### General safety instructions

- The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and must be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids and pharmaceuticals, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids and pharmaceuticals shall be taken into account. • Please observe all further instructions for use of any compatible medical devices used in parallel. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections to users.

### Notes on safe use of product

- The product is designed for use with compatible endoscopes, with generated pressures up to 50 kPa (0.5 bar), by gravity or with a pressure cuff with generated pressures up to 50 kPa (0.5 bar). • Prior

to use, verify the suitability of the product in terms of permissible fluid volume/flow rate with regard to the urethral/uterine anatomy of the patient. • Prior to use, the product shall be checked for damage, tightness, and any misconnections. • Prior to use, existing closure plugs and disconnectable connections must be checked for tightness and retightened if necessary. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70 % or isopropyl alcohol 70 % in combination with chlorhexidine gluconate 2 % can be used for disinfection. Please observe that the exposure time must amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Hang up the container(s) of irrigation fluid at a sufficient height (e.g. approx. 100 cm) above the patient's bladder and ensure there are no kinks to ensure unobstructed fluid inflow by gravity. • A labelling "not for intravenous use" of the prepared irrigation set with irrigation solution is recommended.

### Use of product

Please also observe the exemplary product illustration(s) and all further information/images for product use on and in the product packaging.

General use (cf. exemplary product illustration(s) A/B/C/D):

1. Check if the sterile packaging is damaged or has been unintentionally opened.
2. Remove the product from the sterile packaging and check for possible damage.
3. In case of rigid/semi-rigid containers with irrigation solution, an irrigation set with a vented closure-piercing device (1) (A) is required. In case of bags with irrigation solution, irrigation sets with non-vented (1) (B/C/D) as well as vented closure-piercing devices (1) (A) can be used. Close the ventilation cap of vented closure-piercing device (1) (A) first.
4. Close the roller/closure clamp (6).
5. Remove the protective cap and insert the vented closure-piercing device (1) (A) straight into the uprightly positioned rigid/semi-rigid container with irrigation solution or insert the non-vented closure-piercing device (1) (B/C/D) into the port of the bag with irrigation solution with a rotating movement. Please ensure that the wall of the bag with irrigation solution remains intact during piercing.
6. Fill the drip chamber (4) up to the level ring (3), where applicable (cf. exemplary product illustrations A/B/D):
- 6a. Gently squeeze the drip chamber in case of drip chambers with integrated closure-piercing device (cf. exemplary product illustration A/B).
- 6b. In the case of inline drip chambers, incline the lower part of the drip chamber until the angle between the main line and drip chamber is 45°. The drip chamber can now be filled to the level ring from above (cf. exemplary product illustration D).
7. In the case of rigid/semi-rigid containers with irrigation solution, open the ventilation cap of the vented closure-piercing device (1) (A). For bags with irrigation solution, the ventilation cap remains closed.
8. Open the roller/closure clamp (6) and prime the irrigation set until the terminal connection port (7).
9. Close the roller/closure clamp (6).
10. Connect the connection port (7) of the irrigation set as following:
  - 10a. Step connector (cf. exemplary product illustration A/D (7a)): Insert the step connector into the irrigation port of a transurethral bladder catheter.
  - 10b. Silicone tubing (cf. exemplary product illustration B/D (7b)): Slide the silicone tubing on the irrigation port (nozzle) of the endoscope.
  - 10c. Luer-Lock connector (cf. exemplary product illustration C/D (7c)): Connect the Luer-Lock connector to a Luer-Lock compatible irrigation port of a bladder catheter or an endoscope.
11. Check the connection between the connection port of the irrigation set and the irrigation port of the bladder catheter or endoscope, respectively.
12. Open the roller/closure clamp (6) to start the irrigation.

Use of product configurations with multiple closure-piercing devices (cf. exemplary product illustration D):

The various closure-piercing devices (1) can be used for successive transfer of irrigation solution. Thus, a continuous flow of irrigation fluid

can be ensured. The respective unused line shall be closed with the corresponding closure clamp (2).

For adjustment of the terminal connection ports (7) of the irrigation sets according to the required application, corresponding urology adapters for urological and gynaecological applications can be used, or connection ports (cf. exemplary product illustration B/D) can be disconnected.

To extend the irrigation sets, corresponding extension lines for urological and gynaecological applications can be used.

#### Product compatibility

• The closure-piercing device complies with EN ISO 8536-4 and can be connected with all compatible containers with irrigation solution. • The step connector is compatible with sterile urethral catheters according to EN ISO 20696. • The silicone tubing fits the irrigation port (nozzle) for irrigation fluid of urological/gynaecological endoscopes. • All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/ EN ISO 80369-7.

#### Residual risks

No significant residual risks are known at this time.

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### Period of use

The international and national guidelines for the transfer of fluids and pharmaceuticals shall be taken into account.

A maximum period of use of 24 hours should not be exceeded.

Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

### Product disposal

• The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. • Protective medical gloves must be worn during disposal.

### Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay.

Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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