- en** **Instructions for use - Irrigation set**
for urological applications with cytostatics
- de** **Gebrauchsanweisung - Spülset**
für urologische Anwendungen mit Zytostatika
- cs** **Návod k použití - Iriřační sada**
pro urologická použití s cytostatiky
- da** **Brugsanvisning - Irrigationssæt**
til urologisk anvendelse med cytostatika
- es** **Instrucciones de uso - Set de irrigación**
para aplicaciones urológicas con citostáticos
- fr** **Notice d'utilisation - Set d'irrigation**
pour applications urologiques avec cytostatiques
- it** **Istruzioni per l'uso - Set di irrigazione**
per applicazioni urologiche con citostatici
- nl** **Gebruiksaanwijzing - Irrigatieset**
voor urologische toepassingen met cytostatika
- no** **Brukerveiledning - Irrigasjonssett**
til urologiske behandlinger med cytostatika
- pl** **Instrukcja używania - Zestaw irygacyjny**
do zastosowań urologicznych z użyciem cytostatyków
- pt** **Instruções de utilização - Sistema de irrigação**
para aplicações urológicas com citostáticos
- sv** **Bruksanvisning - Irrigationsaggregat**
för urologiska användningar med cytostatika

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
Irrigation set

for urological applications with cytostatics

Intended purpose

Transfer of fluids and pharmaceuticals for irrigation during urological applications as well as their drainage and body fluids, if applicable

Indications

Therapeutic urological applications in the urinary bladder and prostate that require the transfer of fluids and pharmaceuticals.

Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained staff

Self-application is not permitted and is not considered as intended use.

The current Summary of Product Characteristics of the pharmaceutical manufacturers with regard to incompatibilities (interaction material/ pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg.

The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of urethral anatomy, appropriateness of the solution being used, and the duration of therapy.

User

Qualified medical personnel only. For use in hospitals and specialised practices with conventional environmental and hygiene conditions.

Clinical benefits

The transfer of fluids and pharmaceutical solutions for irrigation in urological applications, as well as their drainage, can bring about clinical benefits that improve the health of the patient.

Irrigation in urological applications carries implicit risks and should therefore only be carried out if it improves the patient's health.

General safety instructions

• The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Sterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and must be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids and pharmaceuticals, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment. • Furthermore, the international and national guidelines for the transfer of fluids and pharmaceuticals shall be taken into account. • Please observe all further instructions for use of any compatible medical devices used in parallel. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • In case of consecutive transfer of different drugs with the same irrigation set, rinsing with a neutral compatible solution is recommended after each medication to avoid potential incompatibility of drug mixtures, maintain the integrity of the line, and prevent microbiologic growth. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections to users.

Notes on safe use of product

• The product is designed for manual application with hand-held

syringes with generated pressures up to 50 kPa (0.5 bar). • Prior to use, verify the suitability of the product in terms of permissible fluid volume/flow rate with regard to the urethral anatomy of the patient. • Prior to use, the product shall be checked for damage, tightness, and any misconnections. • Prior to use, existing closure plugs and disconnectable connections must be checked for tightness and retightened if necessary. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70 % or isopropyl alcohol 70 % in combination with chlorhexidine gluconate 2 % can be used for disinfection. Please observe that the exposure time must amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Excessive forces on Luer-Lock connections can lead to stress cracks.

Use of product

Please also observe the exemplary product illustration(s) and all further information/images for product use on and in the product packaging.

General use (cf. exemplary product illustration(s)):

1. Check if the sterile packaging is damaged or has been unintentionally opened.
2. Remove the product from the sterile packaging and check for possible damage.
3. After checking for possible leakages, remove the prepared connector (min. 2-way, with integrated check valves) connected with the medication/cytostatic syringe(s) by the preparing site out of the transport container. Check if the closure clamps of the connector are closed.
4. Connect a Luer-Lock syringe with 20 ml neutral solution to the connection port for flushing of the connector.
5. Connect the step connector (1) of the irrigation set to the connection port of the bladder catheter as well as the urology adapter (5) to the urine collection bag. Check if the connection is tight.
6. Connect the already prepared connector to the lateral Luer-Lock connection port (3) of the irrigation set.
7. Close the closure clamp (4) of the irrigation set.
8. Open the closure clamp of the medication intended to be used first, as well as the closure clamp for flushing of the connector.
9. Instill the medication from the prefilled medication/cytostatic syringe into the bladder.
10. After complete administration of the medication/cytostatic syringe, flush with the syringe with neutral solution to assure that the Luer-Lock connection port between irrigation set (3) and connector is free of medication/cytostatics.
11. Repeat steps 8 to 10 according to the number of intended medications. Use a new Luer-Lock syringe with 20 ml neutral solution for irrigation every time.
12. Close all closure clamps of the connector.
13. Disconnect the connector from the Luer-Lock connector port (3) of the irrigation set and dispose of according to national as well as internal guidelines. Close the Luer-Lock connection port (3) of the irrigation set with a new closure plug.
14. The medication/cytostatic remains within the bladder depending on the prescribed administration time according to the actual Summary/ Summaries of Product Characteristics of the manufacturer(s).
15. To empty the bladder, open the closure clamp (4) of the irrigation set. The contaminated urine will be collected in the urine collection bag.

Final disconnection of the irrigation set from the bladder catheter in case of single administration/instillation:

16. After the bladder has been emptied completely and in preparation for disconnecting the irrigation set, connect a Luer-Lock syringe with sufficient neutral solution (10-50 ml) to the lateral Luer-Lock connector (3) of the irrigation set. Remove the closure plug beforehand.
17. Close the closure clamp (4) and flush the connection port to the bladder catheter with sufficient neutral solution (5-35 ml according to the respective tubing volume).
18. Close the bladder catheter. The irrigation set with urine collection bag can be disconnected from the bladder catheter and disposed of in accordance with the regulations for the disposal of contaminated medical waste.

19. A new urine collection bag can be connected directly to the bladder catheter for drainage. Open the bladder catheter afterwards.

Final disconnection of the irrigation set from the bladder catheter in case of multiple administrations/instillations and exchanges of the urine collection bag:

16. To replace the urine collection bag, connect a Luer-Lock syringe with sufficient neutral solution (10 ml) to the Luer-Lock connector (3) of the irrigation set. Remove the closure plug beforehand.

17. Close the bladder catheter and flush the connection port to the urine collection bag with approx. 5–10 ml of neutral solution.

18. Then close the closure clamp (4) of the irrigation set. The urine collection bag can be disconnected from the urology adapter (5) and disposed of closed in accordance with the regulations for the disposal of contaminated medical waste.

19. A new urine collection bag can be connected to the urology adapter (5) for drainage. Open the bladder catheter afterwards.

In case of administration of different drugs/cytostatics with the same irrigation set for urological applications, an appropriate manifold connector (multiple lumen) shall be used. After each medication/cytostatic administration, flushing with neutral solution via the connection port for flushing at the manifold connector is recommended.

Product compatibility

• The step connector is compatible with sterile urethral catheters according to EN ISO 20696. • The urological adapter is compatible with

connection ports of urine collection bags according to EN ISO 8669-2.

• All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/EN ISO 80369-7.

Residual risks

No significant residual risks are known at this time.

Period of use

The international and national guidelines for the transfer of fluids and pharmaceuticals shall be taken into account.

A maximum period of use of 24 hours should not be exceeded.

Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

Product disposal

• The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. • Protective medical gloves must be worn during disposal.

Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay.

Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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