



# CODAN

## AeroGuard®

Additional safety for invasive  
blood pressure monitoring

**The decisive connection**



## Blood pressure monitoring sets from CODAN

For many years now, the invasive blood pressure monitoring sets from CODAN have been used in intensive care and anaesthesia to display and register blood pressure continuously and also to enable arterial blood samples to be taken for laboratory medicine. Invasive blood pressure monitoring is performed on far in excess of 50 % of all patients in intensive care in


hospitals and has always been a precise and reliable method for continuously monitoring the cardiovascular system<sup>1,2</sup>.

The name CODAN not only stands for the maintenance of this gold standard, but also for the consistent optimisation of invasive blood pressure monitoring sets in terms of transfer quality, flexibility and patient safety.

## Patient transport

Even during intra-hospital patient transportation, invasive blood pressure monitoring helps medical staff to assess the patient's cardiovascular system. The major challenge faced by medical staff when transporting a patient lies in reliably and continuously maintaining all monitoring and treatment measures for the entire duration of transportation<sup>3</sup>. Consequently, the transportation of intensive care patients who require neurological, respiratory and haemodynamic monitoring is particularly critical and presents a great potential risk, especially in turbulent situations<sup>3</sup>. If a blood pressure monitoring set is stored incorrectly and the roller and sliding clamps of the filling system are not closed, air can, for example, enter the blood pressure monitoring set via the connected saline bag or drip chamber<sup>4</sup>. Air pockets which are guided into a patient's vascular system can cause air embolisms

which, in turn, can have catastrophic consequences for the patient<sup>4,5</sup>. Even just the most minimal amounts of air in the arterial system can not only lead to extremely serious complications, but also major costs for the clinical facility<sup>6</sup>.




### Sources

1. Burchardi H, Larsen R, Marx G, Muhl E, Schölmerich J. Die Intensivmedizin. Dordrecht: Springer; 2011. 1170 p. ISBN: 978-3-642-16928-1. ger. <http://gbv.eblib.com/patron/FullRecord.aspx?p=769894>.
2. Knapp J, Grabowski M, Weidhase L, Bernhard M. Invasive Blutdruckmessung – Schritt für Schritt. *Kardio up*. 2018;14(04):296–303. ger. doi:10.1055/a-0770-3407.
3. Frank O, Bohner K, Dzemali O. Empfehlungen - Innerklinische Transporte kritisch kranker Patienten. Zürich (Schweiz): Patientensicherheit Schweiz, Stiftung für Patientensicherheit; 2014 [accessed 2019 Jul 1]. ger. <https://www.patientensicherheit.ch/>.
4. Holmaas G. Arteriell gassetemboli ved bruk av arteriekran og kontinuerleg intraarteriell blodtrykksmåling. *NAForum*. 2016;29(2):35. nor.
5. Bove A. Arterielle Gasembolie. Lewis Katz School of Medicine, Temple University: MSD Manual; 2017 [accessed 2019 Jul 1]. ger. <https://www.msdmanuals.com/>.
6. Cook LS. Infusion-related air embolism. *J Infus Nurs*. 2013;36(1):26–36. eng. doi:10.1097/NAN.0b013e318279a804.

## AeroGuard® – Providing that extra bit of safety

In order to prevent air from entering into the patient's vascular system at all costs, the safety component AeroGuard® was specially developed for invasive blood pressure monitoring. This filter system is integrated into the classic filling system of the blood pressure monitoring set from CODAN and essentially comprises two chambers and filter membranes each. Due to its hydrophilic features, the 1.2 µm membrane located between the chambers retains all air pockets that would pass through the blood pressure measuring set from the direction of the saline bag. This air is then separated off out of the system via a hydrophobic 0.02 µm membrane.

AeroGuard® therefore not only protects against the undesired penetration of larger amounts of air, but also against tiny air bubbles which could enter the hose system whilst the blood pressure monitoring set is being filled. Moreover, thanks to the 1.2 µm membrane, particles are also filtered out at the same time. Positioning the AeroGuard® above the pressure transducer decouples the filter system from the blood pressure monitoring and, in this way, ensures that the signal transmission is not negatively influenced by the component.




## The outstanding advantages of AeroGuard®

AeroGuard® sets new standards when it comes to patient safety during invasive blood pressure monitoring. The filter system offers:

- Increased patient safety thanks to reliable air separation of the 0.02 µm membrane at ≤ 500 mmHg
- Particle filtration through the 1.2 µm membrane
- Full functionality even after having been wetted externally with NaCl solution or disinfectant\*
- Use of PVC-free materials

\*After completely filling the AeroGuard®

Together with the reliable and precise pressure transducers in the Xtrans® series and the hygienically closed blood sampling systems from CODAN, AeroGuard® offers next-level benefits for all healthcare facilities where invasive blood pressure monitoring sets are used.



# CODAN Worldwide

CODAN is known internationally as a manufacturer and supplier of medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 60 years of research and development. Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

## CODAN Product range

- Infusion sets
- Transfusion sets
- Extension sets & Manifold connectors
- Infusion & Transfusion accessories
- Infusion filters & Filter systems
- Neonatology & Paediatrics products
- Withdrawal, Preparation & Administration systems
- CODAN CYTO®
- Chemoprotect® products
- Urology & Gynaecology products
- Single use syringes
- Invasive Blood Pressure Monitoring
- Infusion Technology Solutions
- Volumetric & Syringe pumps
- Software
- CODAN COMPONENTS

## CODAN Companies

CODAN Medizinische Geräte GmbH · Deutschland  
CODAN pvb Critical Care GmbH · Deutschland  
CODAN pvb Medical GmbH · Deutschland  
CODAN II, S.A. · Portugal  
CODAN US Corporation · USA  
CODAN Inc. · USA  
CODAN NORGE AS · Norge  
CODAN TRIPLUS AB · Sverige  
CODAN Limited · United Kingdom  
CODAN FRANCE Sarl · France  
CODAN Medical AG · Schweiz  
CODAN ARGUS AG · Schweiz  
arcomed AG · Schweiz  
CODAN BV · Nederland  
CODAN s.r.l. · Italia  
CODAN Medical GmbH · Österreich  
CODAN Medical ApS · Danmark  
CODAN DEHA ApS · Danmark  
CODAN MEDITECH s.r.o. · Česká republika  
arcomed Pty Ltd · Australia  
arcomed Ltd · New Zealand

Compliance of the established quality management systems with the provisions of EN ISO 13485, the Council Directive 93/42/EEC and/or Regulation (EU) 2017/745 has been certified by the relevant, competent notified bodies:

## TÜV SÜD Product Service GmbH

CODAN Medizinische Geräte GmbH  
23738 Lensahn, Germany

CODAN pvb Critical Care GmbH  
85661 Forstinning, Germany

CODAN US Corporation  
Santa Ana, CA 92704, USA

CODAN ARGUS AG  
6340 Baar, Switzerland

arcomed AG  
8302 Kloten, Switzerland

## Presafe Denmark A/S

CODAN Medical ApS  
4970 Rødby, Denmark



## Manufacturer

CODAN pvb Critical Care GmbH  
Römerstraße 18  
D-85661 Forstinning  
Tel.: +49 (0) 81 21 · 98 02 0  
codan@codanpvbcc.de  
www.codancompanies.com

**The decisive connection**