



CODAN

Check Valves

The versatile programme of check valves

The decisive connection



CODAN

Check Valves



Small, perfect and safe

Top-class engineering in a mini-format

As with all CODAN products, the CODAN Check Valves also manifest decades of experience in research and development.

Precision in production ensures maximum reliability and ease of movement of the CODAN Check Valve. The valve opens upon the slightest conveying pressure. It closes automatically as soon as a counter-pressure results or the conveying pressure has dropped again.


Even highly viscous infusions pass through the valves securely and reliably.

The CODAN Check Valves have a particularly low priming volume of 0.15 ml and can be integrated into all infusion systems without any problems.

Be it in standard infusion therapy on wards, in the operating theatre, in intensive care or in anaesthetics - CODAN Check Valves can be used in all areas of everyday hospital work.

Integrated safety

CODAN Check Valves are a small detail with maximum importance for patient safety. As a firm or optional integral part in various CODAN products, the check valves reliably prevent reflux or back flow of infusion solution without additional measures and also ensure that blood does not flow back from the patient into the infusion system.



R-LOCK,
Female/Male
Luer-Lock adapter
with CODAN
Check Valve

Inline adapter
with CODAN
Check Valve

Double Y-site
with CODAN
Check Valves

Y-injection site
with CODAN
Check Valve

Double connector
with CODAN
Check Valves

Rotating
Male Luer-Lock
adapter with CODAN
Check Valve

Male Luer-Lock
adapter with CODAN
Check Valve

Versatility with precision

For use in many versions

The matching solution can be found for each application. The CODAN Check Valve can be integrated into an existing system as an individual component. In addition, a variety of connection sets, extensions, connectors, completely pre-fabricated systems and special systems with CODAN Check Valves are available. All variants are available with differing tubing diameters. Connectors with CODAN Check Valves are available as double, triple or quadruple connectors.

Connections

CODAN Check Valves in CODAN infusion set connections make handling particularly safe and simple - the less individual components are used, the more quickly and more reliably the infusion can be prepared and administered. Special importance is given to the Y-shaped injection site in safe and low-risk administration of cytostatics.

Connection variants

The connections provided with a rotating collar which enable secure closure of the Male Luer-Lock without tube torsion are particularly comfortable and user-friendly. Of course CODAN also provides all standard infusion sets with integrated CODAN Check Valves.

CODAN Worldwide

CODAN is known internationally as a manufacturer and supplier of medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 60 years of research and development. Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

CODAN Product range

- Infusion sets
- Transfusion sets
- Extension sets & Manifold connectors
- Infusion & Transfusion accessories
- Infusion filters & Filter systems
- Neonatology & Paediatrics products
- Withdrawal, Preparation & Administration systems
- CODAN CYTO®
- Chemoprotect® products
- Urology & Gynaecology products
- Single use syringes
- Invasive Blood Pressure Monitoring
- Infusion Technology Solutions
- Volumetric & Syringe pumps
- Software
- CODAN COMPONENTS

CODAN Companies

CODAN Medizinische Geräte GmbH · Deutschland
CODAN pvb Critical Care GmbH · Deutschland
CODAN pvb Medical GmbH · Deutschland
CODAN TI, S.A. · Portugal
CODAN US Corporation · USA
CODAN Inc. · USA
CODAN NORGE AS · Norge
CODAN TRIPLUS AB · Sverige
CODAN Limited · United Kingdom
CODAN FRANCE Sarl · France
CODAN Medical AG · Schweiz
CODAN ARGUS AG · Schweiz
arcomed AG · Schweiz
CODAN BV · Nederland
CODAN s.r.l. · Italia
CODAN Medical GmbH · Österreich
CODAN Medical ApS · Danmark
CODAN DEHA ApS · Danmark
CODAN MEDITECH s.r.o. · Česká republika
arcomed Pty Ltd · Australia
arcomed Ltd · New Zealand

Compliance of the established quality management systems with the provisions of EN ISO 13485, the Council Directive 93/42/EEC and/or Regulation (EU) 2017/745 has been certified by the relevant, competent notified bodies:

TÜV SÜD Product Service GmbH

CODAN Medizinische Geräte GmbH
23738 Lensahn, Germany

CODAN pvb Critical Care GmbH
85661 Forstinning, Germany

CODAN US Corporation
Santa Ana, CA 92704, USA

CODAN ARGUS AG
6340 Baar, Switzerland

arcomed AG
8302 Kloten, Switzerland

Presafe Denmark A/S

CODAN Medical ApS
4970 Rødby, Denmark



Manufacturer

CODAN Medizinische Geräte GmbH
Stig Husted-Andersen Straße 11
23738 Lensahn · DE
Tel.: +49 (0) 43 63 · 5111 · Fax: +49 (0) 43 63 · 5112 14
codan@codan.de
www.codancompanies.com

The decisive connection