



CODAN

Connect Z® Clip universal & TRANSPORT-LOCK

Two solutions, one standard: added safety
and reliability in everyday clinical use

The decisive connection




CODAN

Connect Z® Clip universal & TRANSPORT-LOCK

Two solutions, one goal: maximum safety and reliability in the infusion process

CODAN is expanding its portfolio with two practical innovations designed to reliably prevent medication loss and to enhance safety throughout the entire infusion process in a simple and efficient way.

The Connect Z® Clip universal and the TRANSPORT-LOCK increase safety for medical personnel and patients alike – from preparation and transport through to patient applications and safe disposal.



Connect Z® Clip universal Key benefits:

- Secure hold – prevents unintended disconnection during transport, e.g. in pneumatic transport systems or during handling and protects against loss of medication
- One clip for Connect Z® and infusions bags – reduced product variants and simplifies inventory management
- Flexibly adjustable to the infusion bag
- Needle-free mixing adapter remains accessible – even with the fixing clip in place

TRANSPORT-LOCK Key benefits:

- Secure fixation – the infusion set remains reliably connected to the infusion bag
- Increased safety during movement, patient mobilisation or repositioning
- Reduced risk of medication loss
- Reliable support in daily clinical use
- Supports safe and controlled disposal

Connect Z® Clip universal

Universally compatible. Securely connected.

The Connect Z® Clip universal is a universal fixing clip designed to securely connect CODAN Connect Z® secondary lines with infusion bags.

As a further development of our proven Connect Z® Clip variants, the Connect Z® Clip universal simplifies the entire process by requiring only one fixing clip for all compatible infusion bags. This reduces product variants, simplifies inventory management and ensures clear, efficient workflows in pharmacies, practices and hospitals. Detailed information on compatibility is available in the product information and can be provided on request.

By securely fixing the connection, the clip prevents unintended loosening or disconnection – for example during transport, in pneumatic transport systems or during handling. This effectively minimises the risk of accidental medication loss.



TRANSPORT-LOCK

Secure fixation – during transport and throughout infusion

This is a reliable fixing solution designed to securely connect infusion bags with standard infusion sets.

The TRANSPORT-LOCK prevents unintended disconnection of the infusion sets caused by tension or movement – from preparation and transport through the entire infusion process and up to final disposal. As a result, the connection remains securely closed, fluid leakage is avoided and the risk of medication loss is significantly reduced. At the same time, process safety in daily clinical practice is noticeably enhanced.

It can be used with all compatible standard infusion sets. Detailed information on compatibility is available in the product information and can be provided on request.



CODAN Worldwide

CODAN is known internationally as a manufacturer and supplier of medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 60 years of research and development. Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

CODAN Product range

- Infusion sets
- Transfusion sets
- Extension sets & Manifold connectors
- Infusion & Transfusion accessories
- Infusion filters & Filter systems
- Neonatology & Paediatrics products
- Withdrawal, Preparation & Administration systems
- CODAN CYTO®
- Chemoprotect® products
- Urology & Gynaecology products
- Single use syringes
- Invasive Blood Pressure Monitoring
- Infusion Technology Solutions
- Volumetric & Syringe pumps
- Software
- CODAN COMPONENTS

CODAN Companies

CODAN Medizinische Geräte GmbH · Deutschland
CODAN pvb Critical Care GmbH · Deutschland
CODAN pvb Medical GmbH · Deutschland
CODAN TI, S.A. · Portugal
CODAN US Corporation · USA
CODAN Inc. · USA
CODAN NORGE AS · Norge
CODAN TRIPLUS AB · Sverige
CODAN Limited · United Kingdom
CODAN FRANCE Sarl · France
CODAN Medical AG · Schweiz
CODAN ARGUS AG · Schweiz
arcomed AG · Schweiz
CODAN BV · Nederland
CODAN s.r.l. · Italia
CODAN Medical GmbH · Österreich
CODAN Medical ApS · Danmark
CODAN DEHA ApS · Danmark
CODAN MEDITECH s.r.o. · Česká republika
arcomed Pty Ltd · Australia
arcomed Ltd · New Zealand

Compliance of the established quality management systems with the provisions of EN ISO 13485, the Council Directive 93/42/EEC and/or Regulation (EU) 2017/745 has been certified by the relevant, competent notified bodies:

TÜV SÜD Product Service GmbH

CODAN Medizinische Geräte GmbH
23738 Lensahn, Germany

CODAN pvb Critical Care GmbH
85661 Forstinning, Germany

CODAN US Corporation
Santa Ana, CA 92704, USA

CODAN ARGUS AG
6340 Baar, Switzerland
arcomed AG
8302 Kloten, Switzerland

Presafe Denmark A/S

CODAN Medical ApS
4970 Rødby, Denmark



Contact

Global Sales & Marketing Management
CODAN Companies
Tel.: +49 (0) 43 63 · 5111
customerservice@codancompanies.com
www.codancompanies.com

Manufacturer

CODAN Medizinische Geräte GmbH
Stig Husted-Andersen Straße 11
23738 Lensahn · DE
Tel.: +49 (0) 43 63 · 5111 · Fax: +49 (0) 43 63 · 5112 14
codan@codan.de
www.codancompanies.com

The decisive connection