

CODAN

Pain therapy pumps

PCA | PCEA | PIEB | PCS

The decisive connection



NEW
Breath
activated
ArcoAir
available

CODAN

Pain therapy pumps

All-in-One concept

ArcoAir PCA Switch

An ever-increasing number of patients are unable to press a button due to compromised dexterity. arcomed ArcoAir provides a solution to this, with breath-activated PCA. To trigger the patient bolus, the patient simply blows into the mouthpiece and the

programmed patient bolus will be delivered. A bulb-activated attachment for ArcoAir is also available, which can be used to deliver PCA in the MRI Suite (in conjunction with MRIShield).



Chroma μ SP6000 PCA

- Large colour touch screen
- Compact lightweight design, with integrated lock
- Individual PCA protocols
- Paediatric and weight based programming
- Hourly summary of successful and refused demands
- Viewable event logs
- 3 separate lock levels
- Password protected folders





Tailor-made Solutions & a strong Partnership

Creating advanced infusion systems



Made in Switzerland



Innovative Solution



Environment Awareness



Patient Safety

Extensive PCA therapies

Lockable PCA pumps



- Syringe
- Volumetric

PCA Therapies

- Breath-activated PCA, with ArcoAir
- PCA in the MRI, with ArcoAir and MRIShield
- Hyperbaric PCA
- PCS – Patient Controlled Sedation
- PIEB – Programmed Intermittent Epidural Bolus

PCS

Patient Controlled Sedation (PCS) can be used to cover a broad range of areas where sedation is needed, such as burn care, out-patient surgery and endoscopic procedures. Studies show that PCS was preferred by the patients.



CODAN Worldwide

CODAN is known internationally as a manufacturer and supplier of medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 60 years of research and development. Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

CODAN Product range

- Infusion sets
- Transfusion sets
- Extension sets & Manifold connectors
- Infusion & Transfusion accessories
- Infusion filters & Filter systems
- Neonatology & Paediatrics products
- Withdrawal, Preparation & Administration systems
- CODAN CYTO®
- Chemoprotect® products
- Urology & Gynaecology products
- Single use syringes
- Invasive Blood Pressure Monitoring
- Infusion Technology Solutions
- Volumetric & Syringe pumps
- Software
- CODAN COMPONENTS

CODAN Companies

CODAN Medizinische Geräte GmbH · Deutschland
CODAN pvb Critical Care GmbH · Deutschland
CODAN pvb Medical GmbH · Deutschland
CODAN TI, S.A. · Portugal
CODAN US Corporation · USA
CODAN Inc. · USA
CODAN NORGE AS · Norge
CODAN TRIPLUS AB · Sverige
CODAN Limited · United Kingdom
CODAN FRANCE Sarl · France
CODAN Medical AG · Schweiz
CODAN ARGUS AG · Schweiz
arcomed AG · Schweiz
CODAN BV · Nederland
CODAN s.r.l. · Italia
CODAN Medical GmbH · Österreich
CODAN Medical ApS · Danmark
CODAN DEHA ApS · Danmark
CODAN MEDITECH s.r.o. · Česká republika
arcomed Pty Ltd · Australia
arcomed Ltd · New Zealand

Compliance of the established quality management systems with the provisions of EN ISO 13485, the Council Directive 93/42/EEC and/or Regulation (EU) 2017/745 has been certified by the relevant, competent notified bodies:

TÜV SÜD Product Service GmbH

CODAN Medizinische Geräte GmbH
23738 Lensahn, Germany

CODAN pvb Critical Care GmbH
85661 Forstinning, Germany

CODAN US Corporation
Santa Ana, CA 92704, USA

CODAN ARGUS AG
6340 Baar, Switzerland
arcomed AG
8302 Kloten, Switzerland

Presafe Denmark A/S

CODAN Medical ApS
4970 Rødby, Denmark



Contact

Global Sales & Marketing Management
CODAN Companies
Tel.: +49 (0) 43 63 · 5111
customerservice@codancompanies.com
www.codancompanies.com

Manufacturer

arcomed AG
Steinackerstrasse 29
CH-8302 Kloten
info@arcomed.com
www.arcomed.com

The decisive connection